

TECHNICAL DATA SHEET – CALXGLUK


CALXGLUK

5,11,17-tris{2-[4-(β-D-glucopyranosyl)-1H-1,2,3-triazol-1-yl]ethyl}-25-monoheptoxy-26,27,28-trimethoxycalix[4]arene

2018

Information

Compound Name	CALXGLUK	Physical state	Beige powder
Catalogue Number	CALXGLUK_10MG, CALXGLUK_50MG, CALXGLUK_100MG	Purity (HPLC, 280nm)	<95%
Molec. Formula	C ₇₁ H ₉₅ N ₉ O ₂₂	Retention time (RP₁₈ HPLC)^a	t _R = 13.4 min
CAS	nd	CMC	0.025 mM
MW	1426.5 g/mol	Exact Mass	1425.6592
pKa	na		
Percent composition	C, 59.78; H, 6.71; N, 8.84; O, 24.67		
Stability	Store in <-20°C freezer for up to one year		
Solubility	Soluble in water (2.5mM), methanol and DMSO.		
Structure			




References

1- Proprietary patent: "Synthesis of amphiphilic calixarene glycoside detergents and use of same for extracting and stabilizing native functional membrane proteins"; WO 2015158575 A1 ; Publication date: Oct 22, 2015.


TECHNICAL DATA SHEET – CALXGLUK

Biochemical Validation Data



Binding of radioligand on target protein, purified in reference detergent with or without CALXGLUK as an additive

Purified protein was incubated with radioligand in absence (total, blue bars) or presence (Non Specific signal, red bars) of an excess of cold ligand. After filtration on GF/C membranes and washing, scintillation agent was added and radioactivity was detected using a Microbeta2. Specific radioactivity (green bars) corresponds to (total signal) – (non-specific signal).



Stabilization of the Bacteriorhodopsin

After solubilization and purification in OG, reagents were added as stabilizers.

Bacteriorhodopsin's activity was monitored over time at first 25°C and then 37°C by measuring absorbance at 560nm.

Restricted use: Limited Use Label License

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product and progeny, to perform internal research and development for the sole benefit of the purchaser.

Safety: Not known as a hazardous substance or mixture. General industrial hygiene practices must be followed as the use of adapted personal protective equipment for skin and body.

Technical support

For additional product and technical information email our Technical Support team at contact@calixar.com.

Limited product warranty

TO THE EXTENT ALLOWED BY LAW, CALIXAR WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

© 2017 Calixar. All rights reserved. All trademarks are the property of Calixar and its subsidiaries.