

**RECOMBINANT HUMAN ADENOSINE A_{2A} RECEPTOR
 BATCH # AMO_210402**
2021
Protein information

Target Name	Adenosine A _{2A} receptor (A _{2A} receptor)
Catalogue Number	PP1
Class	GPCR Class A
Sequence	Full-length, wildtype sequence, with a N-terminus Strep tag II , 8xHis-tag , and TEV protease cleavage site : MWSHPQFEKHHHHHHHENLYFQG PIMGSSVYITVELAIAVLAILGNVLCWAVWLNSNLQNVNTNY FVVSLLAAADIAVGVLAIPFAITISTGFCAACHGCLFIACFVLVLTQSSIFSLIAIDRYIAIRIPLRYNGL VTGTRAKGIIAICWVLSFAIGLTPMLGWNNCGQPKEGKNHSQGCQVACLFEDVPMNYMVYF NFFACVLVPLLLMLGVYLRIFLAARRQLKQMESQPLPGERARSTLQKEVHAAKSLAIIVGLFALCWLPL HIINCFTFFCPDCSHAPLWLMYLAIVLSHTNSVVPFIYAYRIREFRQTFRKIIRSHVLRQQEPFKAAG TSARVLAAHGSDGEQVSLRLNGHPPGVWANGSAPHPERRPNGYALGLVSGGSAQESQGNTGLPDV ELLSHELKGVCEPPGLDDPLAQDGAGVS
Affinity Tag	His/Strep (both N-terminal)
Origin	Human (Homo sapiens)
Theor. MW	47,7kDa
Accession #	P29274 (UniProt)

Protein production


Expression system	Sf9 insect cells (baculovirus)
Purification	Immobilized Metal Affinity Chromatography
Purity	>90%
Activity	Confirmed by radiobinding assay
Concentration	Up to 5mg/ml
Sample Buffer	25 mM Na ₂ HPO ₄ /NaH ₂ PO ₄ pH 8.0, 150 mM NaCl, 0.05%-0.006% DDM-CHS
Available quantity	From 10µg up to mg scale
References	1- Igonet S et al. Stabilization of native and functional Adenosine receptor. <i>Scientific Reports</i> . 2018 May 25; doi.org/10.1038/s41598-018-26113-0 2- Desuzinges Mandon E. et al. Novel systematic detergent screening method for membrane proteins solubilization. <i>Anal Biochem</i> . 2017 Jan 15;517:40-49.

Miscellaneous

Shipment Temperature	Dry ice
Storage conditions	Store at -80°C




Quality Controls (Purity and Activity):



SDS-PAGE, 4-15% acrylamide gel
Bio-rad Stain-Free™ detection

SDS-PAGE, Stain-Free detection

Purified A_{2A} was migrated on a 4-15% Tris-glycine SDS-PAGE and the total proteins were Stain-Free detected. Black arrows indicate the target. Upper arrows indicate full-length A_{2A}. Lower arrow indicate shorter A_{2A} resulting from partial cleavage of C-term end.




SDS-PAGE, 4-15% acrylamide gel
WB Anti-ICL3 A_{2A} antibody (7F6-G5-A2)

SDS-PAGE, western blotting

Purified A_{2A} was migrated on a 4-15% Tris-glycine SDS-PAGE, transferred to pvdf membrane and immunodetected with a monoclonal anti- ICL3-A_{2A} (7F6-G5-A₂, SCBT). Black arrows indicate the target.

³H-CGS21680 agonist binding to A_{2A}R



QC: Activity measured by radiobinding assay

Binding of [³H]CGS21680 was measured on purified A_{2A}R. A K_D of 194 nM was determined for CGS21680.

Restricted use: All Calixar products or derived products shall be used for internal research only excluding any direct or indirect commercial uses. (Under patents rights – For other uses, please contact us at contact@calixar.com).

Limited Use Label License: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product and progeny, to perform internal research and development for the sole benefit of the purchaser.

Safety: Not known as a hazardous substance or mixture. General industrial hygiene practices must be followed as the use of adapted personal protective equipment for skin and body.

Technical support

For additional product and technical information email our Technical Support team at contact@calixar.com.

Limited product warranty

TO THE EXTENT ALLOWED BY LAW, CALIXAR WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

© 2022 Calixar. All rights reserved. All trademarks are the property of Calixar and its subsidiaries.